

SmartCart by Alphatron

Making the last few metres of the medication round smarter, faster, safer



Introducing The SmartCart

Work smarter, faster and safer on
your medication rounds

At Omnicell we think there is a more effective way for nurses to complete their medication rounds. Working in partnership with Alphatron Medical Innovations, we are transforming your nurses' patient care workflows with the SmartCart, harnessing intelligent technology to deliver:

- Smarter daily medication ordering and management - no note taking against printed lists, instead order and collect medications directly through your Omnicell XT cabinet
- Faster rounds with less repetitive tasks, walking to and from the patient
- Safer delivery of medications with secured individual patient bins ensuring no mixing errors, and transparent closed loop recording of medication administration

By combining the functionality of medication trolleys with computers on wheels, and providing dedicated locking medication bins, the SmartCart transforms medication dispensing rounds for nurses.



Built to last, easy to use



High quality design
and materials



Passive and active
work modes



Reliable fanless
battery system:
runtime 12+ hours



Automatic ergonomic
working position



Electrical height
adjustment for easy
patient interaction



Wifi enabled for access
to EPR, EPMA and
MAR systems



10+ years lifespan

Closing The Loop With Your Omnicell XT Cabinet and The SmartCart

By combining your Omnicell Automatic Dispensing Cabinets (ADC) with the SmartCart you can deliver closed loop medication administration direct to your patient's bedside, whilst realising more value from your existing ADC investments.


Depending on your Electronic Patient Record (EPR) system, there are two ways to integrate your ADCs with the SmartCart. For Epic and Oracle Health users, patient prescription data can be accessed directly from the SmartCart and their medications are queued to the ADC for collection.

For other EPR system users, patient medications can still be ordered and queued to your ADCs with just a few additional clicks.

In both cases, your nursing staff benefit from a single visit to the cabinet during each medication round, whilst patients benefit from extra nurse time at the bedside and the peace of mind that patient safety is your top priority.



Omnicell XT ADC +
AMiS PRO SmartCart



Select medication to dispense
directly to your ADC



SmartCart assigns a medication
bin for each patient, opens bin for
nurse to place all medications



See The SmartCart in Action

To arrange an appointment
email info@omnicell.com

Visit [Omnicell.me](https://omnicell.me) to learn more today.